 The EU directive 93/42/EEC was applied in the design and production of this medical device.

**NAKANISHI INC.** 

[www.nsk-inc.com](http://www.nsk-inc.com)

700 Shimohinata Kanuma-shi  
Tochigi 322-8666,  
Japan

**NSK Europe GmbH** 

[www.nsk-europe.de](http://www.nsk-europe.de)

Elly-Beinhorn-Strasse 8  
65760 Eschborn,  
Germany

**NSK United Kingdom Ltd**

[www.nsk-uk.com](http://www.nsk-uk.com)

Office 5, Gateway1000,  
Arlington Business Park, Whittle Way,  
Stevenage, SG1 2FP, UK

**NSK America Corp**

[www.nskdental.com](http://www.nskdental.com)

1800 Global Parkway  
Hoffman Estates, IL 60192,  
USA

**NSK France SAS**

[www.nsk.fr](http://www.nsk.fr)

32 rue de Lisbonne  
75008 Paris,  
France

**NSK Dental Spain SA**

[www.nsk-spain.es](http://www.nsk-spain.es)

C/ Módena,43 El Soho-Európolis  
28232 Las Rozas, Madrid,  
Spain

**NSK Oceania Pty Ltd**

[www.nskoceania.com.au](http://www.nskoceania.com.au)

Unit 22, 198-222 Young St.  
Waterloo, Sydney,  
NSW 2017, Australia

**NSK Middle East**

[www.nsk-inc.com](http://www.nsk-inc.com)

Room 6EA-701, 7th Floor, East Wing No.6  
Dubai Airport Free Zone,  
PO Box 54316 Dubai, UAE

**NSK Asia Pte Ltd**

[www.nsk-inc.com](http://www.nsk-inc.com)

1 Maritime Square,  
#09-33 HarbourFront Centre,  
Singapore 099253

Specifications are subject to change without notice.

'11.09.00 

OPERATION MANUAL

**NSK**

 0197

Please read this Operation Manual carefully  
before use, and file for future reference.

High Speed Air Turbine Handpiece

S-Max 

S-Max MT

Powerful Partners®

**MADE IN JAPAN**

OM-T0345E 003

<b>EN</b>	Use Environment	Temperature	0 – 40 °C	Transportation and Storage environment	Temperature	-10 – 50 °C
		Humidity	30 - 75% HR		Humidity	10 - 85% RH
		Atmospheric Pressure	700 - 1,060hPa		Atmospheric Pressure	500 - 1,060hPa

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**Intended Use** : This handpiece is designed only for dental treatment use.

## CAUTION

- When operating the handpiece always consider the safety of the patient.
- Do not attempt to disassemble the handpiece nor tamper with the mechanism.
- Check for vibration, noise and overheating outside the patient's mouth. If any abnormalities are found do not operate. Contact authorized dealer for service.
- Should the handpiece not function normally, cease operation immediately and return the handpiece to your authorised NSK Dealer for service.
- Depressing the chuck push button system while the bur is rotating could result in OVERHEATING of the handpiece head. Caution must be exercised during use to keep cheek tissue away from the chuck push button system. Contact with soft tissue may cause the push button to depress and burn injury to the patient may occur.
- Avoid impact on the handpiece. Do not drop the handpiece.
- Remove the bur only after the handpiece has completely stopped rotating.
- Do not use long surgical burs. Do not use burs longer than 26mm.
- This product is classified as Class 1 LED Product. (M500WLED, M500BLED, M600WLED, M600BLED)
- Do not fixate the eyes of LED light.(M500WLED, M500BLED, M600WLED, M600BLED)
- If any abnormalities (dark, do not light on or flashing) are found for LED, stop using the product immediately and contact authorized dealer. (M500WLED, M500BLED, M600WLED, M600BLED)
- Use the power source with the supply capability of 15W or less for LED lighting.
- Use a power source which meets the following requirements. (M500WLED, M500BLED, M600WLED, M600BLED)
  1. The electricity supply of the power source is below 15W both under normal and single failure conditions.

2. The power source uses a SELV circuit for electricity supply.
3. The output voltage of the power source is within the range recommended by the manufacturer of this product.

## 1. Specification

Model	M500L	M500	M500KL/M500SL/ M500WL/M500BL	M500K	M500WLED M500BLED	M600L	M600	M600KL/M600SL/ M600WL/M600BL	M600K	M600WLED M600BLED	
Head Type	Miniature Head					Standard Head					
Rotation Speed	380,000~450,000min <sup>-1</sup>					360,000~430,000min <sup>-1</sup>					
Drive Air Pressure	0.20MPa (2.0kgf/cm <sup>2</sup> ) ~0.25MPa (2.5kgf/cm <sup>2</sup> )		0.25MPa (2.5kgf/cm <sup>2</sup> ) ~0.3MPa (3.0kgf/cm <sup>2</sup> )			0.20MPa (2.0kgf/cm <sup>2</sup> ) ~0.25MPa (2.5kgf/cm <sup>2</sup> )		0.25MPa (2.5kgf/cm <sup>2</sup> ) ~0.3MPa (3.0kgf/cm <sup>2</sup> )			
Chuck Type	Push Button Chuck										
Burs	Size: ISO1797-1 Type3 ø1.59~1.60mm										
	Friction Grip Short/Standard Shank					Friction Grip Standard Shank					
Optic	Glass Rod	-	Glass Rod	-	White LED	Glass Rod	-	Glass Rod	-	White LED	
Voltage	-					AC/DC 3.3±0.05V		-			AC/DC 3.3±0.05V
Activate Voltage	-					2.8 - 4.0 V		-			2.8 - 4.0 V
Consumption current	-					0.38A Typical (3.3V)		-			0.38A Typical (3.3V)

## 2. Connection / Disconnection to the Coupling

First make sure the manufacturer's coupling is firmly attached to the handpiece tubing and the drive air pressure is correct.

### CONNECTION (Fig.1)

To connect the handpiece to the coupling align the handpiece to the coupling and push together (per the coupling manufacturers instructions).

### DISCONNECTION (Fig.1)

Release the handpiece from the coupling (per the coupling manufacturers instructions) and remove the handpiece.

**NSK COUPLING** : Pull back Connector Ring to release the handpiece.


Fig. 1

### 3. Insertion & Removal of the bur

#### TO INSERT THE BUR (Fig. 2)

Depress the Push Button and insert the bur into the chuck until it is secure. Release the button.

#### TO REMOVE THE BUR (Fig. 3)

Stop the handpiece. Depress the Push Button firmly and remove the bur.


Fig. 2


Fig. 3

#### CAUTION

Test that the bur is secure by gently pulling and pushing the bur without depressing the Push Button.

### 4. Information about burs

- Do not use non-standard burs. The ISO standard shank diameter is  $\varnothing 1.59 - \varnothing 1.60$  mm.
- Do not use bent, worn, damaged, or non-concentric burs. Such burs can cause damage to the handpiece.
- Do not exceed the bur speed recommended by the bur manufacturer.
- Always keep the bur shank clean. Entry of hard debris into the chuck via the bur shank could cause rotation slip and also prevent the bur from being securely located in the chuck.
- Fully depress the Push Button and insert the bur into the chuck until it is secure.
- Do not use short shank burs in Standard Head handpieces.

### 5. Lubrication, Cleaning and Sterilization

Refer to the Maintenance Guide **OM-Z0135E**.

All handpieces can be autoclavable up to Max.135°C.


All handpieces can be cleaned and disinfected with a Thermo-Disinfector.

## 6. Replacing the optic illumination bulb (M500WL, M600WL, M500BL, M600BL Series)

Remove the handpiece from the coupling. Loosen and remove the Taper Ring. (Fig. 4). Pull out the joint then remove the bulb holder, then the bulb. Align the pins of the new bulb with the points in the bulb holder and insert the bulb. Place the bulb holder back into the joint. (Fig. 5).

Firmly tighten the Taper Ring.

*Optional Bulb : TA Bulb (Pack of 3) Order No.Y900132*


\* Figure shows M500WL/,M600WL

Fig. 4

### CAUTION

- Do not touch the glass section of new bulb.
- Screw the taper ring tightly. A loose fitting taper ring could result in water and air leakage. It is easier to screw in the taper ring when the bulb holder is lightly pressed in by hand.

<M500WL, M600WL>

<M500BL, M600BL>


Fig. 5

## 7. Replacing the O-rings (M500L, M500, M600L, M600 Series)

Replace the O-rings if water is present in the exhaust air line. This is an indication of possible water leakage within the coupling. ALWAYS change the complete set of O-rings. Remove the handpiece from the coupling. Loosen and remove the Taper Ring at the rear of the handpiece. (Fig. 4) Gently remove each O-ring with your thumb. Insert the complete set of new O-rings in the correct grooves. (Fig. 6). Replace and firmly tighten the Taper Ring.

*Optional O-ring : PTL O-ring Set (Pack of 5) Order No.Y900580*


Fig. 6

### CAUTION

- Make certain that the Taper Ring is firmly tightened. If the Taper Ring is loose water and air leakage could occur.
- Do not disassemble the taper-ring except the O-ring changing.

## 8. Cleaning the cellular glass optic rod (M500L, M600L, M500KL, M600KL, M500SL, M600SL, M500WL, M600WL, M500WLED, M600WLED, M500BL, M600BL, M500BLED, M600BLED)

Wipe clean the cellular glass optic rod entry point (Fig. 7) and exit point (Fig. 8) with an alcohol-immersed cotton swab. Remove all debris and oil.


Fig. 7


Fig. 8

### CAUTION

Do not use a sharp tool to clean the cellular glass optic rod. It could damage the glass and reduce the light transmission. If illumination becomes dim please contact your dealer.

## 9. Replacing the Turbine Cartridge

- (1) Remove the turbine cartridge from the handpiece
  - 1) Insert a bur into the chuck.
  - 2) Locate the head cap wrench tool on the head cap then turn the wrench counter-clockwise. Remove the cap.
  - 3) Use the bur to gently lever the entire turbine cartridge out of the head.
  - 4) Remove debris from inside the empty head.
- (2) Insert a new NSK turbine cartridge
  - 1) Gently insert the new NSK turbine cartridge into the head section, always making certain that the turbine cartridge locating pin is aligned in the slot inside the head.
  - 2) Use fingers to screw the head cap back into place and then use the cap wrench to further tighten until the cap is securely in place.

*Optional Cartridge : SX-MU03 (Cartridge for Miniature head) Order No.P1001250*  
*SX-SU03 (Cartridge for Standard head) Order No.P1002250*


Fig. 9

## CAUTION

- Use only the genuine NSK-manufactured turbine cartridge. Do not use repaired or re-built NSK turbine cartridge. Non-conformance to this shall immediately void the manufacturer's warranty, and the manufacturer shall have no liability to the user and the patient for any incident that may arise from use of such turbine cartridge.
- The head cap screw thread is very fine. To prevent damage to the thread always first tighten the handpiece head cap using fingers, then secure with the head cap wrench.
- Make sure that the two O-rings (Fig.10) are correctly located on the cartridge.
- Use only NSK genuine cartridge. If you use our cartridge, specification may not be satisfied and we can not warranty of this situation.


Fig. 10

## 10. Disposing product

Consult with dealer from whom you purchased it about waste disposal.

\*KaVo® and MULTiflex® are registered trademarks of Kaltenbach & Voigt GmbH & Co., Germany. Sirona® is a registered trademark of Sirona Dental Systems GmbH, Germany.

W&H® and Roto Quick® are registered trademarks of W&H Dentalwerk Bürmoos GmbH, Austria. Bien-Air® and Unifix® are registered trademarks of Bien-Air Dental S.A., Switzerland.

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